



DATA ANALYSIS AND VALIDATION ON ONCOKB, GATEWAYSEQ, AND MYELOSEQ-HD

Justin Caringal & Ajay Khanna, Dr. Eric Duncavage Lab
Washington University of Saint Louis
Department of Pathology & Immunology



ANALYSIS OF ONCOKB API CALLS IN GATEWAYSEQ

STARTING POINT AND GOALS

- GatewaySeq: A tumor-only, high coverage targeted next generation sequencing assay for the identification of gene mutations, copy number alterations, microsatellite instability, tumor mutational burden, and gene fusions
- Evaluate three different methods for looking up variants in OncoKB from MSKCC via web API for interpretation provided to physicians

- byGenomicChange, byProteinChange, byHGVSg

BRAF →

7,140453136,140453136,A,T

p.V600E

7:g.140453136A>T

<https://pathologyservices.wustl.edu/items/gatewayseq-ngs-panel-with-interpretation/>




TABLE ANNOTATIONS [142 SAMPLES]


Comparisons	Differences
<i>Genomic vs. Protein</i>	+60
<i>Genomic vs. HGVSg</i>	-2
<i>Protein vs. HGVSg</i>	-60

	Successful Hits	Failed Hits
<i>Genomic</i>	538	1214
<i>Protein</i>	484	1268
<i>HGVSg</i>	540	1212
<i>Total Searches</i>	1752	


byGenomicChange Calls



byProteinChange Calls



byHGVSg Calls



DISCREPANCIES BETWEEN GATEWAYSEQ AND ONCOKB

- AKT1
- ATM
- B2M
- CD79B
- CHEK1
- DICER1
- FGFR1
- HRAS
- ITPKB
- MYD88
- NF1
- NFKBIE
- NTRK3
- PAX8
- PTPRD
- RAD51B
- RSPO3
- SGK1
- SMARCA4
- SMARCB1
- TCF3
- TFEB
- TNFAIP3

Gene/Transcript	1752
<i>Correct ID</i>	1489
<i>Incorrect ID</i>	263
<i>Discrepancy Frequency</i>	0.15

- OncoKB and GatewaySeq differed in choice of gene transcript
- GatewaySeq uses the Ensembl Canonical transcript
- AA coordinates may differ, could affect lookups

GATEWAYSEQ CONCLUSION

- Overall, small differences between API calls
- byProteinChange lagged behind both byGenomicChange and byHGVSg
- In cases of successful hit discrepancies
 - byProteinChange failure: No p.syntax available (e.g. splice variant)
 - byGenomicChange failure: Complex variants may not be found
- HGVSg is marginally better than Genomic
 - Example: ARID1A, ENST00000324856 (PASS) TAG→AA, complex variant
- Future Directions: Might be better to use byHGVSg in the pipeline



COVERAGE COMPARISONS IN MYELOSEQ-HD

MYELOSEQ-HD

- Targeted sequencing assay for 49 genes and gene hotspots that are recurrently mutated in myeloid neoplasms, such as MDS and AML
- Uses a high coverage UMIs-based error corrected sequencing approach to achieve >95% sensitivity for previously identified mutations with VAFs $\geq 0.25\%$

<https://pathologyservices.wustl.edu/items/myeloseq/>



READ COVERAGE

- New MyeloSeq-HD feature:
 - For previously identified variants, limit of detection (LOD) depends on sampling error
 - Sampling error depends on coverage
- Paired-End Sequencing reads can overlap with small enough fragment sizes
 - What is the coverage of loci in overlaps, 1 (collapsed) or 2?

>>>>>>>>>>>>>>>>>>>

-----X-----

read 1

fragment

read2

x: locus



GENERAL PIPELINE

- MyeloSeq-HD coverage:
 - Dragen aligns FASTQ files, paired-end sequencing (150 base pairs), reports region coverage → BED
 - Call specific variants (Dragen, Pindel, combination) → VCF/JSON
 - Pindel & Dragen-Pindel uses custom Python script
- Goal: Compare output VCF and BED coverage to measure discrepancies



DRAGEN [67 VARIANTS]

JSON-BED

2496

2496

486

-9 to 0 (57)

-10

-11

-20

-78

-1664

-1778

-3865

- Coverage: order of thousands
- Most (57/67) were near-identical (-9 to 0)
- Examined NRAS variant alignments

CHROM	POS	REF	ALT	GENE	JSON	BED	JSON-BED
chr1	114716127	C	T	NRAS	5664	7442	-1778

- With Samtools mpileup (-Q10 -q20), able to reproduce JSON (variant caller) coverage of 5664
- May want to generate own coverage, or look at options for Dragen BED (region) coverage



JSON-BED

4629
2683
1734
1683
1624
650
615
578
473
445
432
286
174
98
11
0 to 5 (11)
-33
-3056
-4187

DRAGEN-PINDEL (29 VARIANTS)

- Many more large discrepancies
 - Tends to be higher in JSON (variant caller)
 - Possibly due to filtering script
- Did not have time to investigate further
 - We are thinking about removing filtering script

Pindel-only contains 2 cases





THANK YOU!

Special thanks to the Spencer Lab!

Questions?

CREDITS: This presentation template was created by
Slidesgo, including icons by **Flaticon**, and infographics &
images by **Freepik**



